

Customized production is available.

Free standing cabinet is designed for enough storage space.
Construction with heavy duty stainless steel 304

SPECIFICATION

All SUS304 1.0t & 1.2t plate
Glass : 5t strong glass
Shelf : adjustable 30mm intervals
Adjustable foot
Key locking
Requested size available

AVAILABLE OPTIONS


1. Hinged or Sliding door
2. Drawer
3. Sliding work table
4. Bin storage
5. ABS basket
6. Casters



INSTRUMENT CABINET

Model	Size (L×D×H / mm)	Remark
KS-IC3	900×450×1800	

Hinged glass door with key locking
4 adjustable shelves
High adjustable foot



INSTRUMENT CABINET

Model	Size (L×D×H / mm)	Remark
KS-IC4	900×450×1800	Open type
KS-IC4-1	1200×450×1800	Open type

4 adjustable shelves
High adjustable foot



INSTRUMENT CABINET

Model	Size (L×D×H / mm)	Remark
KS-IC1	1200×500×1800	

Sliding glass doors with key locking
Sliding work table
Plastic 24 bins
2 adjustable shelves
High adjustable foot



INSTRUMENT CABINET

Model	Size (L×D×H / mm)	Remark
KS-IC2	1200×500×1800	

Sliding glass doors with key locking
Sliding work table
3 SUS drawers
3 adjustable shelves
High adjustable foot



INSTRUMENT CABINET

Model	Size (L×D×H / mm)	Remark
KS-IC5	900×450×1800	

Hinged glass door with key locking
Hinged SUS door
2 SUS drawers
2 adjustable shelves
High adjustable foot



INSTRUMENT CABINET

Model	Size (L×D×H / mm)	Remark
KS-IC6	1200×500×1800	

Sliding glass doors with key locking
Sliding SUS door
3 adjustable shelves
High adjustable foot



INSTRUMENT CABINET

Model	Size (L×D×H / mm)	Remark
KS-IC15	1300×550×1800	

Sliding glass doors with key locking
Sliding work table
ABS BASKET
SUS wire basket
2 adjustable shelves
High adjustable foot



INSTRUMENT CABINET

Model	Size (L×D×H / mm)	Remark
KS-IC16	1300×550×1800	

Sliding glass doors with key locking
Sliding work table
ABS BASKET
SUS wire basket
2 adjustable shelves
High adjustable foot



INSTRUMENT CABINET

Model	Size (L×D×H / mm)	Remark
KS-IC19	900×450×1800	

Sliding work table
Plastic drawers
2 SUS drawers
3 adjustable shelves
High adjustable foot



INSTRUMENT CABINET

Model	Size (L×D×H / mm)	Remark
KS-IC20	1300×550×1800	


Sliding glass doors with key locking
Sliding work table
Plastic 24 bins
ABS BASKET
SUS wire basket
High adjustable foot



INSTRUMENT CABINET

Model	Size (L×D×H / mm)	Remark
KS-IC17	1300×550×1800	

Sliding glass doors with key locking
ABS BASKET
SUS wire basket
High adjustable foot



INSTRUMENT CABINET

Model	Size (L×D×H / mm)	Remark
KS-IC18	700×550×1800	

ABS BASKET
SUS wire basket
High adjustable foot



INSTRUMENT CABINET

Model	Size (L×D×H / mm)	Remark
KS-IC21	1200×500×1800	




Sliding glass doors with key locking
Sliding work table
4 adjustable shelves
High adjustable foot



INSTRUMENT CABINET

Model	Size (L×D×H / mm)	Remark
KS-IC22	1200×500×1800	

Sliding glass doors with key locking
4 adjustable shelves
High adjustable foot



►KS-PMS-100

►KS-PMS-200

►KS-PMS-300

SPECIAL FEATURE



1. Attach UV Lamp check LED Lamp to the front of the equipment



2. Splash bump design



3. Install pollution control soap case

APPLICATION

- Surgical operating room
- Delivery room
- Emergency room
- ICU room
- Radiology room
- Other places in need of sterilized water

FEATURE

- Heavy duty stainless steel 304 construction for rust proof and corrosion resistance
- Sterilization system by UV lamps and dual filterings by main filter 5um and pre-filter 10um
- Sterilized water of over 99% is provided
- Easy operation by knee-water and soap
- Water pressure control device

BASIC OPERATION AND FUNCTION

- 1.Heating
 - 1) Hot water supply in winter
 - 2) Fast warming by one heater (2KW)
 - 3) Two safety thermo-switches to control the in-creasin pressure of the warming tank
 - 4) Constant water temperature by safety thermo-switches
- 2.Sterilization
 - 1) The water can be sterilized to 99.5%(0.6T/H) by UV water sterilizing lamps
 - 2) The UV lamps are made of dual quart tubes for electrical safety
 - 3) Lamp us time can be read as life-timer is attached
- 3.Filtering
 - 1) pre-filter 10um
 - 2) Main-filter 5um
 - 3) Thorough filtering by dual-filters
 - 4)Easy filter replacement

SCRUB STATION

Model	Size (L×D×H / mm)	Soap Pump
KS-PMS-100	1000×600×1150	1(220V-3W)
KS-PMS-200	1400×600×1150	1(220V-3W)
KS-PMS-300	2000×600×1150	2(220V-3W)

- Treatment : Germs removal by filtration
- Heater : 2.0Kw×1EA
- Water : 20~60℃
- Water time : Adjustable Within 30sec
- Filter : Pre-Filter 1EA, Main-Filter 1EA
- Sterilizer : Ultra Violelet Ray Flow-Sterilizing Systme (15W x 2EA)
- U.V Lamp : 2,000 Hours, Using Time indicates on control panel with digital (Filter life time check)
- Water supply : Min 1.0kg/cm² Max 5.0kg/cm²
- Water temperature control : Thermostatic Mixing Valve
- Soap filling check : Microprocessor Auto Control
- Power source : AC 220V, 50/60Hz
- Power consumption : 2.2Kw
- Water tank capacity : Approx 30L
- Net weight [kg] 75 / 105 / 135



SINK CANINET

Model	Size (L×D×H / mm)	Remark
KS-SK01	1200×600×(850+100)	-
KS-SK02	1500×600×(850+100)	-
KS-SK03	1800×600×(850+100)	-

- SUS 304 seamless welded construction
- Hinged SUS door
- 3 drawers
- 1 faucet
- 1 sink bowl
- 100mm back plate on the top
- Height adjustable
- Customized size or spec. available



WALL-MOUNTED INSTRUMENT CABINET

Model	Size (L×D×H / mm)	Remark
KS-WMC01	3700×500×2300	Top Cabinet

■ As well as this model type, we can make a customized product and install them according to the request of the customer.

COMPOSITION UNITS

- ① Basket Cabinet Accessories
 - ABS Basket (Low Type) 2EA
 - ABS Basket (High Type) 4EA
 - SUS Wirt Basket (Low Type) 2EA
 - SUS Wirt Basket (High Type) 2EA
 - Multi Box 24EA (Partition)
 - Supporting Slide Table.
- ② Instrument Cabinet Accessories.
 - Shelf 2EA
- ③ Instrument Cabinet Accessories.
 - Shelf 2EA
 - Multi Plastic Bin 24EA(Partition)
- ④ Supporting Slide Table
- ⑤ The Top Cabinet



WALL-MOUNTED INSTRUMENT CABINET

Model	Size (L×D×H / mm)	Remark
KS-WMC03	3400×500×2300	Top Cabinet



WALL-MOUNTED INSTRUMENT CABINET

Model	Size (L×D×H / mm)	Remark
KS-WMC04	3200×500×2300	Top Cabinet



WALL-MOUNTED INSTRUMENT CABINET

Model	Size (L×D×H / mm)	Remark
KS-WMC02	4500×500×2300	Top Cabinet

■ As well as this model type, we can make a customized product and install them according to the request of the customer.


COMPOSITION UNITS

- ①② Instrument Drawer Cabinet
 - Shelf 2EA
- ③ Instrument Bin Cabinet Accessories.
 - Multi Box (Low) 34EA (Partition)
 - Multi Box (High) 32EA (Partition)
- ④ Instrument Cabinet Accessories.
 - Shelf 2EA
- ⑥ The Top Cabinet
- ⑦ Drawer (Swivels Sliding Rail)



WALL-MOUNTED INSTRUMENT CABINET

Model	Size (L×D×H / mm)	Remark
KS-WMC05	4200×500×1800	-



WALL-MOUNTED INSTRUMENT CABINET

Model	Size (L×D×H / mm)	Remark
KS-WMC06	2400×500×2300	Top Cabinet

■ As well as this model type, we can make a customized product and install them according to the request of the customer.

COMPOSITION UNITS

- ① Basket Cabinet Accessories
 - ABS Basket (Low Type) 2EA
 - ABS Basket (High Type) 4EA
 - SUS Wirt Basket (Low Type) 2EA
 - SUS Wirt Basket (High Type) 2EA
 - Multi Box 24EA (Partition)
 - Supporting Slide Table.
- ② Instrument Cabinet Accessories.
 - Shelf 2EA
- ③ The Top Cabinet



RACK BASKET SYSTEM

Model	Size (L×D×H / mm)	Remark
KS-RBS	2450×550×1800	STS Wine Basket 560×410×125(H)



PACS CABINET


Model	Size (L×D×H / mm)	Remark
KS-PC01	1200×500×2300	Top Cabinet

■ Option : Monitor, Viewbox (2Panel)



RACK SYSTEM

Model	Size (L×D×H / mm)	Remark
KS-RACK100	consultation	



INSTRUMENT TABLE (ROUND TYPE)

Model	Size (L×D×H / mm)	Remark
KS-ITR1	1200×550×900	Round
KS-ITR2	1600×550×900	Round
KS-ITR3	1900×550×900	Round

- SUS 304 construction
- Round shaped 2 shelves with 3 side guard rail
- Round rubber bumper against wall
- Six 4 inch casters (two with brakes)
- Customized size or spec. available



INSTRUMENT TABLE (I-TYPE)

Model	Size (L×D×H / mm)	Remark
KS-ITI1	900×500×900	flat type
KS-ITI2	1200×500×900	flat type
KS-ITI3	1600×500×900	flat type

- SUS 304 construction
- Bottom shelf with 3 side guard rail
- Round rubber bumper against wall
- Four 4 inch casters (two with brakes)
- Customized size or spec. available



INSTRUMENT TABLE

Model	Size (L×D×H / mm)	Remark
KS-ITC3S	625×450×1000	—



INSTRUMENT TABLE (V-TYPE)

Model	Size (L×D×H / mm)	Remark
KS-ITV1	1200×500×900	Type V
KS-ITV2	1600×500×900	Type v
KS-ITV3	1900×500×900	Type v

- SUS 304 construction
- V-shaped 2 shelves with 3 side guard rail
- Round rubber bumper against wall
- Six 4 inch casters (two with brakes)
- Customized size or spec. available



INSTRUMENT TABLE (G-TYPE)

Model	Size (L×D×H / mm)	Remark
KS-ITG1	900×500×900	guard
KS-ITG2	1200×500×900	guard
KS-ITG3	1600×500×900	guard

- All Stainless Steel 304 Construction
- 4 " Swivel Brake Caster 4EA
- Two Self shaped with guard rail



INSTRUMENT TABLE

Model	Size (L×D×H / mm)	Remark
KS-ITC2D	625×425×940	—



MAYO STAND

Model	Size (L×D×H / mm)
KS-MSF	490×330×850~1350
KS-MSF-1	600×400×850~1350

- Operating room use



MAYO STAND

Model	Size (L×D×H / mm)
KS-MSH	490×330×850~1350
KS-MSH-1	600×400×850~1350

- Outpatient Clinic use



MAYO STAND(AUTO)

Model	Size (L×D×H / mm)
KS-MSA	490×330×850~1350
KS-MSA-1	600×400×850~1350

-



WASH BASIN STAND

Model	Size (L×D×H / mm)	Remark
KS-WBS01	350×350×900	Single

- All Stainless Steel Pipe Frame
- Stainless Basin
- 2 " Caster 4EA



WASH BASIN STAND

Model	Size (L×D×H / mm)	Remark
KS-WBS02	700×350×900	Double

- All Stainless Steel Pipe Frame
- Stainless Basin
- 2 " Caster 4EA



MAYO STAND

Model	Size (L×D×H / mm)
KS-02-043	490×330×850~1350

-



MAYO STAND (HAND WORK)

Model	Size (L×D×H / mm)
KS-MSFA	490×330×850~1350
KS-MSFA-1	600×400×850~1350




MAYO TRAY

Model	Size (L×D×H / mm)
KS-MSA-01-B	600×400×18




MAYO TRAY2

Model	Size (L×D×H / mm)
KS-MSA-B	490×330×18



RACK DRYING GLOVE

Model	Size (L×D×H / mm)	Remark
KS-RDG20	600×400×1500	20쪽
KS-RDG30	800×400×1500	30쪽
KS-RDG50	1200×400×1500	50쪽



SUS WIRE BASKET

Model	Size (L×D×H / mm)	Remark
KS-SWB-1	560×410×130	SUS Bar
KS-SWB-2	560×410×190	SUS Bar



SUCTION BOTTLE 3000 CC

Model	Size (L×D×H / mm)	Remark
KS-SB-30C	φ350×500	Cap include

- Plastic 3000cc bottle (Polycarbonate)
- Sus bottle cap : 2 air line
- Option accessories : Silicon horse
- Overflow function



SUCTION BOTTLE HOLDER

Model	Size (L×D×H / mm)	Remark
KS-SB-H80C	8000cc	Regularator

- Swivel caster 2 " 4EA
- Urethane bumper
- Plastic 8000cc bottle (Polycarbonate)
- Sus bottle cap : 2 air line
- Silicon horse (2M-2EA)



SUCTION BOTTLE HOLDER(Dual)

Model	Size (L×D×H / mm)	Remark
KS-SB-H30CD	3000cc	Regularator



SUCTION BOTTLE 8000 CC

Model	Size (L×D×H / mm)	Remark
KS-SB-80C	φ350×550	Bottle
KS-SB-80CA	—	Cap

- Plastic 8000cc bottle (Polycarbonate)
- Sus bottle cap : 2 air line
- Option accessories : Silicon horse
- Overflow function



SUCTION BOTTLE HOLDER

Model	Size (L×D×H / mm)	Remark
KS-SB-H30C	3000cc	Regularator

- Swivel caster 2 " 4EA
- Urethane bumper
- Plastic 3000cc bottle (Polycarbonate)
- Sus bottle cap : 2 air line
- Silicon horse (2M-2EA)



CHEST BOTTLE STAND

Model	Size (L×D×H / mm)	Remark
KS-LBS-02	165×145×285	—

- All sus304 construction



KICK BUCKET

Model	Size (L×D×H / mm)	Remark
KS-KB02	φ350×350	—

- SUS 304 construction
- Round rubber bumper
- Stainless steel bucket
- Four 2 inch swivel casters



KICK BUCKET

Model	Size (L×D×H / mm)	Remark
KS-KB05	φ430×900	—

- SUS 304 construction
- Round rubber bumper
- Stainless steel bucket
- Four 2 inch swivel casters



CHEST BOTTLE STAND

Model	Size (L×D×H / mm)	Remark
KS-LBS-01	180×100×285	—

- All sus304 construction



KICK BUCKET

Model	Size (L×D×H / mm)	Remark
KS-KB03	φ350×700	Handle Type


- SUS 304 construction
- Round rubber bumper
- Handle
- Stainless steel bucket
- Four 2 inch swivel casters



KICK BUCKET

Model	Size (L×D×H / mm)	Remark
KS-KB04	φ350×140	—

- SUS 304 construction
- Round rubber bumper
- Four 2 inch swivel casters



TUBE STAND

Model	Size (L×D×H / mm)	Remark
KS-TS1	210×170×500	STS304



TUBE STAND

Model	Size (L×D×H / mm)	Remark
KS-TS2	300×200×180	STS/AL



SIDE LAMP

Model	Size (L×D×H / mm)	Remark
KS-SL88A	500×1250×1850	45,000 LUX



INSTRUMENT CART

Model	Size (L×D×H / mm)	Remark
KS-MO4-D6	1160×500×1110	



INSTRUMENT CART

Model	Size (L×D×H / mm)	Remark
KS-MO5-D9	1700×500×1110	



MONITOR CART

Model	Size (L×D×H / mm)	Remark
KS-MO1	600×600×1300	—



MONITOR CART

Model	Size (L×D×H / mm)	Remark
KS-MO2	550×500×1300	—



MONITOR CART

Model	Size (L×D×H / mm)	Remark
KS-MO3	600×500×1200	—

■ Possibility to make a customized cart.



INSTRUMENT CART

Model	Size (L×D×H / mm)	Remark
KS-ITC-4W	635×475×1170	



INSTRUMENT CART

Model	Size (L×D×H / mm)	Remark
KS-ITC-4D	635×475×1170	