


[3 MOTOR]

KS-DM-3



HOSPITAL ELECTRIC MOTION BED

Model	Size (L×D×H / mm)	Remark
KS-DM-3	2180×940×400~700	3 motor

제품 컬러



리모컨




기능

PP 블로우사출 매트리스지지판

PP 블로우 사출형매트리스판이 견고하고 위생적으로 사용이 가능하여 반영구적 사용이 가능



도면



제품정보

침대 높이	최저 400mm 최고 700mm (매트리스 제외)
침대 길이	최대 2180mm
침대 폭	940mm
안전작동하중	200kg
매트리스	820×1910×80T
사이드레일	알루미늄 사이드레일
매트리스 지지판	PP사출
식탁	멤퍼오일 가스 쇼바 장착

[2 MOTOR]

KS-DM-2



HOSPITAL ELECTRIC MOTION BED

Model	Size (L×D×H / mm)	Remark
KS-DM-2	2180×940×500~600	2 motor

제품 컬러



리모컨



선택사양



유로하우징 5인치 캐비바퀴

선택사양



낙상방지 접이식 식탁 칸막이


기능

PP 블로우사출 매트리스지지판

PP 블로우 사출형매트리스판이 견고하고 위생적으로 사용이 가능하여 반영구적 사용이 가능



도면




제품정보

침대 높이	최저 500mm 최고 600mm (매트리스 포함)
침대 길이	최대 2180mm
침대 폭	940mm
안전작동하중	180kg
매트리스	820×1910×80T
사이드레일	알루미늄 사이드레일
매트리스 지지판	PP사출
식탁	멤퍼오일 가스 쇼바 장착

[1 MOTOR]

KS-DM-1



HOSPITAL ELECTRIC MOTION BED

Model	Size (L×D×H / mm)	Remark
KS-DM-1	2180×940×500~600	1 motor

제품 컬러



리모컨



선택사양



유로화우정 5년저 커버배워

선택사양



낙상방지 잠이식 식탁 칸막이


기능

PP 블로우사출 매트리스지지판

PP 블로우 사출형매트리스판이 견고하고 위생적으로 사용이 가능하여 반영구적사용이 가능



도면



제품정보

침대 높이	최저 500mm 최고 600mm (매트리스 포함)
침대 길이	최대 2180mm
침대 폭	940mm
안전작동하중	180kg
매트리스	820×1910×80T
사이드레일	알루미늄 사이드레일
매트리스 지지판	PP사출
식탁	멤퍼오일 가스 쇼바 장착

[2 CRANK]

KS-DM-2A



HOSPITAL ELECTRIC MOTION BED

Model	Size (L×D×H / mm)	Remark
KS-DM-2A	2180×940×500~600	2 crank

제품 컬러




기능

핸들로 등판 및 다리판 각도 조절 가능



도면




제품정보

침대 높이	최저 500mm 최고 600mm (매트리스 포함)
침대 길이	최대 2160mm
침대 폭	940mm
안전작동하중	180kg
매트리스	820×1910×80T
사이드레일	알루미늄 사이드레일
매트리스 지지판	철판타공

[1 CRANK]

KS-DM-1A



HOSPITAL ELECTRIC MOTION BED

Model	Size (L×D×H / mm)	Remark
KS-DM-1A	2180×940×500~600	1 crank

제품 컬러



그린

스카이블루


오크

기능

핸들로 등판 및 다리판 각도 조절 가능
PP 블루시울 매트리스 지지판




도면



제품정보

침대 높이	최저 500mm 최고 600mm (매트리스 포함)
침대 길이	최대 2180mm
침대 폭	940mm
안전작동하중	180kg
매트리스	820×1910×80T
사이드레일	알루미늄 사이드레일
매트리스 지지판	PP블루시울



Operating Room
Appliances

Nursing
Appliances

Consulting Room
Appliances

Supply Room
Appliances

Hospital Bed

Carts

Monitor Carts

IV Pole Stand

Monitor Arm

Wall Mount

Pole Mount

Accessory



PILLOW

Model	Size (L×D×H / mm)	Remark
KS-PL-01	further consultation	–

Requested size available



MATRESS

Model	Size (L×D×H / mm)	Remark
KS-MAT01	1800×700×50	General (50T)
KS-MAT02	1870×880(830)×50	for children (50T)

Requested size available



EXAMINING TABLE

Model	Size (L×D×H / mm)	Remark
KS-KY-026	1800×600(700)×650	–


Powder coated steel frame
40T sponge with artificial leather cover
Leg : 30 x 50 mm steel pipe
Foldable aluminum side rail at both side



EXAMINING TABLE

Model	Size (L×D×H / mm)	Remark
KS-KY-027	1800×600(700)×650	–

Powder coated steel frame
40T sponge with artificial leather cover
Leg : 30 x 50 mm steel pipe
Foldable aluminum side rail at both side
5 inch casters (Two with brakes)



ATTENDANT BED

Model	Size (L×D×H / mm)	Remark
KS-KY-014	1600×550×250	–

Powder coated steel frame
Two 2 inch casters



EXAMINING TABLE

Model	Size (L×D×H / mm)	Remark
KS-KY-025	1800×600(700)×650	–

Powder coated steel frame
40T sponge with artificial leather cover
Leg : 30 x 50 mm steel pipe



EXAMINING BED

Model	Size (L×D×H / mm)	Remark
KS-EB	1800×700×650	Steel

Powder coated steel frame
40T sponge with artificial leather cover
Leg : 30 x 50 mm steel pipe



EXAMINING BED

Model	Size (L×D×H / mm)	Remark
KS-EBD	1800×700×650	Drawer Type

Powder coated steel frame
40T sponge with artificial leather cover
Leg : 30 x 50 mm steel pipe
Three wood drawers



CABINET (WOOD)

Model	Size (L×D×H / mm)	Remark
KS-KY-056	400×400×700	1 Person



EXAMINING BED

Model	Size (L×D×H / mm)	Remark
KS-EB300	1900×600×650	Steel

80T sponge with artificial leather cover



SHOWER TROLLEY (STANDARD)

Model	Size (L×D×H / mm)	Remark
KS-STC2	1900×600×750	Standard

SPECIFICATION

Stainless Steel Pipe 304 Construction.
Top plate : Long hole plate.
Side rail : swing Type 2EA.
4Inch Swivel Caster 4EA, (Brake 2EA)



When I lowered the guide rail



CABINET (WOOD)

Model	Size (L×D×H / mm)	Remark
KS-KY-057	450×420×1800	1 Person



CABINET (WOOD)

Model	Size (L×D×H / mm)	Remark
KS-KY-058	450×420×1800	Personal Rocker



INFANT BASSINET STAND

Model	Size (L×D×H / mm)	Remark
KS-IBS1	720×440×800	High Quality

- Stainless Steel Pipe 304 Frame
- P.V.C Clear crib with 50mm Mattress
- 4 Inch Swivel Caster 4EA



INFANT BASSINET STAND

Model	Size (L×D×H / mm)	Remark
KS-IBS2	720×440×800	Standard

- Stainless Steel Pipe 304 Frame
- P.V.C Clear crib with 50mm Mattress
- 4 Inch Swivel Caster 4EA



STRETCHER CART

Model	Size (L×D×H / mm)	Remark
KS-SC1	1900×600×700	Standard

- Stainless Steel 304 Construction
- 5 " Swivel Brake Caster 4EA
- Accessories : Mattress, I.V Pole



STRETCHER CART

Model	Size (L×D×H / mm)	Remark
KS-SC2	1900×600×700	Isolated passive

- Stainless Steel 304 Construction
- 5 " Swivel Brake Caster 4EA
- Accessories : Mattress, I.V Pole



SHOWER TROLLER

Model	Size (L×D×H / mm)	Remark
KS-STC	2000×750×700	—

FEATURE & USAGE

1. It can be easily bathed for patients or people who are in poor motion.
2. Convenient SIDE ONE-TOUCH type SIDERAIL.
3. Minimize corrosion due to all STAINLESS STEEL production
4. Convenient to use by attaching drains to a tarpaulin

FUNCTIONING (Side Rail)



3단 STRETCHER CART

Model	Size (L×D×H / mm)	Remark
KS-SC3	1900×700×750	Steel

- Body : Steel (Powder Painting Management)
- 5 " Swivel Brake Caster 4EA
- Accessories : Mattress, I.V Pole



3단 STRETCHER CART

Model	Size (L×D×H / mm)	Remark
KS-SC4	1900×600×750	SUS

- Stainless Steel 304 Construction
- 5 " Swivel Brake Caster 4EA
- Accessories : Mattress, I.V Pole



STRETCHER CART (1단/토탈락/유압높낮이)

Model	Size (L×D×H / mm)	Remark
KS-HDSC200	1825×615×540~870	—


- Stainless Steel 304 Construction
- 5 " Swivel Brake Caster 4EA
- Accessories : Mattress, I.V Pole



STRETCHER CART(MRI)

Model	Size (L×D×H / mm)	Remark
KS-HDSC200	1825×615×540~870	—

- Stainless Steel 304 Construction
- 5 " Swivel Brake Caster 4EA
- Accessories : Mattress, I.V Pole



3단 PROTECTOR BED / CHAIR

Model	Size (L×D×H / mm)	Remark
KS-YBA1003	561×1777×338	BED SIZE
	561×865×750	CHAIR SIZE



3단 PROTECTOR BED / CHAIR

Model	Size (L×D×H / mm)	Remark
KS-YBA1004	500×1800×340	BED SIZE
	500×865×750	CHAIR SIZE



STRETCHER CART (2단/개별락/유압높낮이)

Model	Size (L×D×H / mm)	Remark
KS-HDSC200	1825×615×540~870	Weight 55 kg



Central locking caster type (Option)