



IV POLE STAND

Model	Size (LxDxH / mm)	Remark
KS-IV-BHR-U	φ 550 × 1250~2000	—

SUS 304 pipe
Outer pipe 27 pie
Inner pipe 22 pie
4 hooks (SUS 304)
ABS five leg base
2.5 inch single wheel
Adjustabl height by knob
ABS plastic handle
Mesh basket
bottle holder



IV POLE STAND

Model	Size (LxDxH / mm)	Remark
KS-IV-BHR-UC	φ 550 × 1250~2000	—

SUS 304 pipe
Outer pipe 27 pie
Inner pipe 22 pie
4 hooks (SUS 304)
ABS five leg base
2.5 inch single wheel
Adjustabl height by knob
Circle plastic basket with handle
Bottle holder



IV POLE STAND

Model	Size (LxDxH / mm)	Remark
KS-IV-BHR	φ 550 × 1250~2000	—

SUS 304 pipe
Outer pipe 27 pie
Inner pipe 22 pie
4 hooks (SUS 304)
ABS five leg base
2.5 inch single wheel
Adjustabl height by knob
ABS plastic handle
Mesh basket



IV POLE STAND

Model	Size (LxDxH / mm)	Remark
KS-IV-HR	φ 550 × 1250~2000	—

SUS 304 pipe
Outer pipe 27 pie
Inner pipe 22 pie
4 hooks (SUS 304)
ABS five leg base
2.5 inch single wheel
Adjustabl height by knob
ABS plastic handle



IV POLE STAND

Model	Size (LxDxH / mm)	Remark
KS-IV-BR-U	φ 550 × 1250~2000	—

SUS 304 pipe
Outer pipe 27 pie
Inner pipe 22 pie
4 hooks (SUS 304)
ABS five leg base
2.5 inch single wheel
Adjustabl height by knob
ABS plastic handle
bottle holder



IV POLE STAND

Model	Size (LxDxH / mm)	Remark
KS-IV-NR-U	φ 550 × 1250~2000	—

SUS 304 pipe
Outer pipe 27 pie
Inner pipe 22 pie
4 hooks (SUS 304)
ABS five leg base
2.5 inch single wheel
Adjustabl height by knob
bottle holder



IV POLE STAND

Model	Size (LxDxH / mm)	Remark
KS-IV-BR	φ 550 × 1250~2000	—

SUS 304 pipe
Outer pipe 27 pie
Inner pipe 22 pie
4 hooks (SUS 304)
ABS five leg base
2.5 inch single wheel
Adjustabl height by knob
Mesh basket



IV POLE STAND

Model	Size (LxDxH / mm)	Remark
KS-IV-NR-C	φ 550 × 1250~2000	—

SUS 304 pipe
Outer pipe 27 pie
Inner pipe 22 pie
4 hooks (SUS 304)
ABS five leg base
2.5 inch single wheel
Adjustabl height by knob
Circle plastic basket with handle



IV POLE STAND

Model	Size (LxDxH / mm)	Remark
KS-IV-HR-U	φ 550 × 1250~2000	—

SUS 304 pipe
Outer pipe 27 pie
Inner pipe 22 pie
4 hooks (SUS 304)
ABS five leg base
2.5 inch single wheel
Adjustabl height by knob
Mesh basket
bottle holder



IV POLE STAND

Model	Size (LxDxH / mm)	Remark
KS-IV-NR	φ 550 × 1250~2000	—

SUS 304 pipe
Outer pipe 27 pie
Inner pipe 22 pie
4 hooks (SUS 304)
ABS five leg base
2.5 inch single wheel
Adjustabl height by knob



IV POLE STAND

Model	Size (LxDxH / mm)	Remark
KS-IV-N	φ 550 × 1270	—

SUS 304 pipe
Outer pipe 27 pie
Inner pipe 22 pie
ABS five leg base
2.5 inch single wheel
Adjustabl height by knob



Foldable hooks (Option)



Five leg base



Height adjustable knob



Circle ABS basket with handle



IV POLE STAND

Model	Size (L×D×H / mm)	Remark
KS-IV-A-BHR-U	φ 550 × 1250~2000	—

SUS 304 pipe
Outer pipe 27 pie
Inner pipe 22 pie
4 hooks (SUS 304)
AL five leg base
2.5 inch single wheel
Adjustabl height by knob
ABS plastic handle
Mesh basket
bottle holder



IV POLE STAND

Model	Size (L×D×H / mm)	Remark
KS-IV-A-BHR-UC	φ 550 × 1250~2000	—

SUS 304 pipe
Outer pipe 27 pie
Inner pipe 22 pie
4 hooks (SUS 304)
AL five leg base
2.5 inch single wheel
Adjustabl height by knob
Circle plastic basket with handle
Bottle holder



IV POLE STAND

Model	Size (L×D×H / mm)	Remark
KS-IV-A-BHR	φ 550 × 1250~2000	—

SUS 304 pipe
Outer pipe 27 pie
Inner pipe 22 pie
4 hooks (SUS 304)
AL five leg base
2.5 inch single wheel
Adjustabl height by knob
ABS plastic handle
Mesh basket



IV POLE STAND

Model	Size (L×D×H / mm)	Remark
KS-IV-A-HR	φ 550 × 1250~2000	—

SUS 304 pipe
Outer pipe 27 pie
Inner pipe 22 pie
4 hooks (SUS 304)
AL five leg base
2.5 inch single wheel
Adjustabl height by knob
ABS plastic handle



IV POLE STAND

Model	Size (L×D×H / mm)	Remark
KS-IV-A-BR-U	φ 550 × 1250~2000	—

SUS 304 pipe
Outer pipe 27 pie
Inner pipe 22 pie
4 hooks (SUS 304)
AL five leg base
2.5 inch single wheel
Adjustabl height by knob
ABS plastic handle
bottle holder



IV POLE STAND

Model	Size (L×D×H / mm)	Remark
KS-IV-A-NR-U	φ 550 × 1250~2000	—

SUS 304 pipe
Outer pipe 27 pie
Inner pipe 22 pie
4 hooks (SUS 304)
AL five leg base
2.5 inch single wheel
Adjustabl height by knob
bottle holder



Model	Size (L×D×H / mm)	Remark
KS-IV-A-BR	φ 550 × 1250~2000	—

SUS 304 pipe
Outer pipe 27 pie
Inner pipe 22 pie
4 hooks (SUS 304)
AL five leg base
2.5 inch single wheel
Adjustabl height by knob
Mesh basket



IV POLE STAND

Model	Size (L×D×H / mm)	Remark
KS-IV-A-NR-C	φ 550 × 1250~2000	—

SUS 304 pipe
Outer pipe 27 pie
Inner pipe 22 pie
4 hooks (SUS 304)
AL five leg base
2.5 inch single wheel
Adjustabl height by knob
Circle plastic basket with handle



IV POLE STAND

Model	Size (L×D×H / mm)	Remark
KS-IV-A-HR-U	φ 550 × 1250~2000	—

SUS 304 pipe
Outer pipe 27 pie
Inner pipe 22 pie
4 hooks (SUS 304)
AL five leg base
2.5 inch single wheel
Adjustabl height by knob
Mesh basket
bottle holder



IV POLE STAND

Model	Size (L×D×H / mm)	Remark
KS-IV-A-NR	φ 550 × 1250~2000	—

SUS 304 pipe
Outer pipe 27 pie
Inner pipe 22 pie
4 hooks (SUS 304)
AL five leg base
2.5 inch single wheel
Adjustabl height by knob



IV POLE STAND

Model	Size (L×D×H / mm)	Remark
KS-IV-A-N	φ 550 × 1270	—

SUS 304 pipe
Outer pipe 27 pie
Inner pipe 22 pie
AL five leg base
2.5 inch single wheel
Adjustabl height by knob



Foldable hooks (Option)



Five leg base



Height adjustable knob



Circle ABS basket with handle